



Monoart

Euronda Monoart®

EN

The original!



Euronda®

Monoart

Kit Monoart®

100%

Monoart® Disposable Dental Kit offers complete clinical protection with premium Italian quality and coordinated 10-colour design. Manufactured by Euronda, the kit combines safety, hygiene, and functionality for professional dental use.

Each kit includes:

- 100 Towel Up! towels
- 100 tissues
- 100 saliva ejectors
- 100 plastic cups
- 100 Protection 3 face masks
- 80 PG30 waterproof aprons



Colours

Blue, blue lagoon, green, light blue, lilac, lime, orange, pink, white, yellow



Size

Box 40 x 25 x 19 cm (lxdxh)



Weight

3,2 kg



Packaging

1 box containing 1 kit



Towel Up!

Monoart® Towel Up! disposable towels are manufactured in Italy using multi-layer pure cellulose tissue bonded with a waterproof polyethylene backing for maximum absorbency and protection. Designed without adhesives, the towels ensure improved hygiene and durability in clinical environments. The resealable 100% recyclable packaging protects against dust and contamination while allowing single-towel dispensing to reduce cross-contact.

Available accessories:

- White towel dispenser (16.5 × 16.5 × 12.5 cm)
- Adjustable brass and chrome-plated steel towel chain (approx. 45 cm)



Colours

Black, blue, blue lagoon, burgundy, grey, light blue, light green, lilac, lime, natura, orange, pink, white, yellow



Size

Pack 17 x 8 x 11 cm (lxdxh)
Towel 33 x 45 cm



Packaging

Pack of 50 towels with 'open/close' system
Carton of 10 packs (500 towels in total)

Practical

Compact and practical design for improved clinic organization and workflow efficiency. Innovative packaging ensures hygiene, easy access, and availability in all Monoart® colour variants.

Safe

Manufactured without adhesives using sustainably sourced pure cellulose and non-toxic polyethylene for safe, hygienic, and eco-conscious clinical use.

Hygienic

Produced on automated manufacturing lines and sealed in protective packaging to maintain maximum hygiene and minimize contamination risk during storage and clinical use.

PG30 aprons

Disposable patient protection aprons supplied in roll format, manufactured from multi-layer absorbent paper with embossed waterproof polyethylene backing, produced without adhesives for improved hygiene and material integrity. Designed for clinical comfort, the aprons feature enlarged coverage and ergonomic neck and shoulder shaping for enhanced patient protection. Pre-printed dashed alignment zones support accurate and efficient placement during use. Made in Italy.

Product options:

- Monoart® PG30 apron mix: 6 rolls in 6 colours (blue, lilac, lime, orange, pink, yellow)

Available accessory:

- Roll apron dispenser container (34.3 × 13.2 × 14.2 cm)



Colours

Blue, blue lagoon, burgundy, light blue, light green, lilac, lime, natura, orange, pink, white, yellow



Size

Roll of 80 aprons:
apron 61 x 53 cm (lxd)
Roll of 60 aprons:
apron 81 x 53 cm (lxd)



Packaging

Pack of 6 rolls

Protective

Designed with ergonomic contouring, the aprons provide full coverage, including shoulders, and stay securely in place around the neck for complete patient protection.

Practical

Equipped with a tie system for secure neck fastening and perforated sections for easy, single-piece use.

Absorbent

Glue-free, mechanically embossed aprons ensure strong layer adhesion and faster liquid absorption for enhanced patient protection.

PG20 aprons

Patient protection aprons on a roll, made from waterproof polyethylene with perforated sections for easy single-use. Made in Italy.



Waterproof

Constructed from polyethylene film for full waterproofing and reliable patient protection.

Protective

Specially contoured aprons provide full-body coverage, including shoulders, and stay securely stretched when tied around the neck.

Practical

The apron features a neck tie system and perforated sections for easy, single-piece use.



Colours

Blue, green, white



Size

Roll of 200 aprons:
apron 81 x 53 cm (lxd)



Packaging

Pack of 6 rolls

Baby towels

Roll of highly absorbent, multi-ply cellulose towels with waterproof polyethylene backing, designed for smaller patients.

Available accessory:

- Adjustable brass and chrome-plated towel chain (approx. 45 cm).



Functional

Highly absorbent, multi-ply cellulose towels with embossed waterproof polyethylene backing for effective protection.

Practical

Compact roll format provides a space-efficient and practical alternative to traditional folded towels for clinical use.

Suitable for small children

Available in two bright colours, designed in a reduced size to ensure effective protection and comfort for pediatric patients.



Colours

Light blue, pink



Size

Roll of 100 towels:
towel 35 x 40 cm



Packaging

Pack of 24 rolls

Tissues

Soft and durable, made from 100% double-ply cellulose (31 g/m²) sourced from FSC-certified forests managed under strict environmental, social, and economic standards. Made in Italy.



Practical

Practical dispenser packaging ensures quick access and easy single-sheet dispensing for efficient clinical use.

Resistant

Constructed from pure double-ply cellulose fibre wadding, offering high strength and superior absorbency for clinical use.

Delicate

Monoart® tissues are soft-touch and skin-friendly, designed for comfortable clinical use.



Colour

White



Size

Pack 23 x 12 x 5 cm (lxdxh)
Towel 20 x 21 cm



Packaging

Single pack (100 tissues)
Box of 40 packs

Saliva ejectors

Flexible disposable saliva ejectors made from non-toxic transparent PVC with an integrated wire core to maintain bending angle stability and prevent leakage. Equipped with a smooth rounded fixed or removable tip for patient comfort, especially in cases of sensitive mucosa. Fully manufactured in Italy under strict quality control using automated, low-contact production systems.

Product mix:

- Monoart® saliva ejector assortment: 12 bags (100 pcs each) in 6 colours (blue, lilac, lime, orange, pink, yellow)

Accessory:

- Autoclavable saliva ejector adaptor (up to 134°C), compatible with 1.6 cm connections.

Eco-friendly

Saliva ejectors are produced on water-recirculating lines to reduce waste, with eco-friendly packaging: recyclable plastic bags and boxes made from certified paper.

Safe

Made from non-toxic, phthalate-free PVC using fully automated production, with multiple quality checks to ensure high standards and minimal defects.

For every need

Offered in two lengths, with fixed or removable rounded tips designed to protect sensitive mucous membranes.



Classification

Medical Device Class IIa



Colours

Blue, blue lagoon, burgundy, green, lilac, lime, orange, pink, transparent, yellow



Size

6.5 mm x 12,5 cm (Øxl)
6.5 mm x 15 cm (Øxl)



Variants

With fixed or removable cap



Packaging

Bag of 100 pieces
Box of 50 bags
(total 5000 pieces)



CE Ref

Monoart® Aspirator
EM15, EM15 Saliva Ejector,
EM15Sp eichelsauger

Aspirator tips

Monoart® aspirator tips, made in Italy, come in three sizes: EM19 (11.4 cm), EM21 (12.9 cm), and EM40 (20 cm). EM19 and EM21 feature ergonomic non-slip grips, rounded tips, and compact design for patient comfort and quiet, hygienic use. EM40 is ideal for hard-to-reach areas.

Accessory:

- EM40 autoclavable tip adaptor (up to 134°C) to convert 1.6 cm connections to 1.1 cm outlets.



Classification

Medical Device Class IIa



Colours

EM19 & EM21: blue, burgundy, green, grey, lilac, lime, orange, pink, yellow
EM40: grey



Size

EM19: 1.6 cm x 11,4 cm (Øxl)
EM21: 1.6 cm x 12,9 cm (Øxl)
EM40: 1.1 cm x 20 cm (Øxl)



Packaging

Pack of 10 pieces
in 1 colour



CE Ref

Monoart® EM19
Monoart® EM21
Monoart® Cannula EM40,
EM40 Aspirator Tip, EM40
Absaugkanüle

Rubber dam

Medium-thickness latex dental dam (0.18–0.23 mm) offers enhanced tear resistance and natural antibacterial properties. Ideal for endodontic and restorative procedures, it isolates the treatment area and reduces aerosol spread, supporting infection control.



Insulating

Provides a reliable barrier between the treatment site and the mouth, minimizing bacterial contamination.

Effective

Soft and flexible, it ensures a dry, stable working area while enhancing visibility and procedural efficiency.

Functional

When combined with double surgical suction, it minimizes aerosol production, reducing the risk of viral and bacterial transmission.



Classification

Medical Device Class I



Colours

Green, grey



Size

15x15 cm



Packaging

Pack of 36 pieces



CE Ref

Monoart Latex Dental Dam

Cotton rolls

Pure long-fibre cotton rolls with high absorbency (up to five times their weight) and low bioburden (<140 CFU/pc) for safe clinical use.

Accessory:

- Cotton roll dispenser, available in white and light blue.



Classification

Medical Device Class I



Size

0.8 cm x 3.7 cm (Øxl)
1 cm x 3.7 cm (Øxl)
1.2 cm x 3.7 cm (Øxl)



Packaging

0.8 cm x 3.7 cm:

Pack of 300 g
Box of 30 packs

1 cm x 3.7 cm:

Pack of 300 g
Box of 60 packs

1.2 cm x 3.7 cm:

Pack of 300 g
Box of 30 packs



CE Ref

Monoart®RULLIDI
COTONE, COTTON
ROLLS, WATTEROLLEN

Functional

Absorbent and flexible, these cotton rolls maintain their shape and position while providing comfort and durability.

Versatile

Offered in three diameters (0.8, 1, and 1.2 cm) to suit different patients and clinical or surgical needs

Eco-friendly

Provides a reliable barrier between the treatment site and the mouth, minimizing bacterial contamination.

Plastic cups

Disposable dental cups made of non-toxic polypropylene with rounded rims and 200 cc capacity, embossed with the Euronda logo for quality assurance. Available in multiple colours and four floral patterns for clinic personalization.

Product assortments:

- Mix Monoart® 1000: 1000 cups in 6 colours
- Mix Monoart® 3000: 3000 cups in 6 colours

Accessory:

- Wall-mounted plastic cup dispenser, holds up to 70 coloured 200 cc cups or 80 white 166 cc cups.

Ergonomics

Monoart® plastic cups feature rounded rims and an ergonomic grip for enhanced patient comfort.

Coloured

Available in 12 colours and 4 floral patterns, adding style and vibrancy to the dental clinic.

Safe

Constructed from certified non-toxic polypropylene to ensure hygiene and safety.



Colours

Black, blue, blue lagoon, floral light blue, floral lilac, floral lime, floral pink, fuchsia, green, light blue, lilac, lime, orange, pink, red, white, yellow



Capacity

200 cc (all colour)
166 cc (white only)



Packaging

Pack of 100 pieces
Box with 30 packs

Paper cups

Foglie

Monoart® Foglie paper cups are made from sustainably sourced European cellulose and produced in Europe to reduce CO₂ emissions. They provide the same water resistance and safety as plastic cups, offering an eco-friendly alternative.



Colour

White with leaf pattern in shades of green and ecru



Capacity

200 cc



Packaging

Pack of 50 pieces
Box with 20 packs

Functional

Monoart® Foglie paper cups offer an eco-friendly alternative to plastic cups while maintaining the same quality, performance, and reliability.

Original

A fresh design adds a personal touch to the Natura line, complementing Monoart® colours: lime, green, white, and ecru.

Eco-friendly

Made entirely from paper sourced from forests managed in accordance with strict environmental, social and economic standards.

Sanitizer gel

Hydroalcoholic hand gel formulated for frequent use without rinsing. Contains 70% alcohol by weight (76% by volume) for rapid reduction of germs and bacteria while helping keep hands soft, moisturized, and pleasantly scented. Made in Italy.



Effective

Designed for dental clinics, it quickly cleans and disinfects, providing effective protection against germs and bacteria.

Functional

Convenient dispenser design; water-free, residue-free, and non-sticky for easy use.

Delicate

Dermatologically tested, lightly scented, and gentle, leaving skin soft, hydrated, and suitable for daily use.



Classification

Cosmetic product
CPNP reference:
1862270



Format

500 cc



Ingredients

Denatured alcohol, Water, Glycerine, Carbomer, Parfum, Aminomethyl propanol, Geraniol, Linalool, Limonene



Packaging

Carton of 6 pieces

Instant ice

Monoart® instant ice pack features an advanced formulation for prolonged and consistent cooling performance. Its compact size ensures precise localized application, while the cellulose-based pouch enhances tactile sensitivity and patient comfort. Manufactured in Italy using automated, hands-free processes to ensure high quality and safety standards.



Classification

Medical Device
Class IIa



Size

15 x 17.5 cm



Shelf life

5 years



Packaging

Pack of 24 pouches



CE Ref

Monoart®Ghiaccio Istantaneo, Instant Ice, Kältekomresse

Gentle on the skin

Designed for sensitive areas, it conforms to the treatment site and feels gentle on the skin for optimal patient comfort.

Compact

Compact pouch design delivers targeted cooling directly to the treatment area.

Durable

Pressure-activated nitrogen salt mixture ensures prolonged and high-quality cooling.

Various accessories



Towel dispenser

Practical white dispenser for disposable towels. Dimensions: 19.5x13.5x17 cm (lxdxh)



Towel chain

Adjustable chrome-plated brass towel chain with steel crocodile clips. Length: 45 cm



Apron dispenser

Practical dispenser container for apron roll. Dimensions: 34x13x14 cm (lxdxh)



Saliva ejector adaptor

Class IIa Medical Device autoclavable to 134°C. Can be used to adapt from a connection with a diameter of 1.6 cm to an outlet with a diameter of 6.5 mm. Pack of 10 pieces.

CER ef: Monoart® Adattatore per Aspirasaliva, Adapter for Saliva Ejector, Adapter für Speichelsauger



Aspirator tip adaptor

Class IIa Medical Device autoclavable to 134°C. Can be used to adapt from a connection with a diameter of 1.6 cm to an outlet with a diameter of 1.1 cm. Pack of 10 pieces.

CE Ref: Monoart® Adattatore per Cannula, Adapter for Aspirator Tip, Adapter für Absaugkanüle



Cotton roll dispenser

Practical dispenser for cotton rolls, available in white and light blue. Dimensions: 11.5x5.5x10 cm (lxdxh)



Cup dispenser

Practical plastic cup dispenser, wall-hanging, holds up to 70 coloured Monoart® cups (200 cc) and up to 80 white cups (166 cc). Dimensions: 11 x 32 cm (Øxh)

Floral line

Add colour and style to your clinic with Euronda Monoart®. The Floral line introduces a fresh spring design in four colours (plus one for masks), featuring polypropylene cups and four-ply floral face masks. Combines safety, protection, and distinctive Monoart® style.



Floral cups

Floral cups – 200 cc disposable cups made of non-toxic polypropylene, durable and reliable, available in four vibrant floral colours.



Face masks Protection Floral

Protection Floral face masks – four-ply protective masks with a floral pattern in five colours, featuring soft rounded elastics and an extended double-track nose bridge for comfort and fit.



Classification

Face masks Protection Floral: Medical Device Class I
EN 14683: type IIR
ASTM F 2100: level 3



Colours

Floral cups: floral light blue, floral lilac, floral lime, floral pink

Face mask Protection Floral: floral burgundy, floral light blue, floral lilac, floral lime, floral pink



CE Ref

Monoart® Face Mask Protection 3 Floral

Face masks Protection 3

Disposable medical face masks designed to prevent airborne microorganism transmission. High breathability (ΔP 25.4 Pa/cm²) exceeds EN 14683 standards, with a particulate and antibacterial filter (99.7% efficiency) and built-in nose bridge.

Three-layer construction:

- Outer: soft, coloured, water-repellent Spunbond fabric
- Middle: Meltblown non-woven filter
- Inner: hydrophilic Spunbond for sweat and moisture absorption

Available with soft side elastics (all colours) or adjustable ties (light blue).

Product mix:

- Mix Monoart® Protection 3 – 6 packs (50 pcs each) in blue, lilac, lime, orange, pink, and yellow.



Classification

Medical Device Class I
EN 14683: type II
ASTM F 2100: level 1



Colours

Blue lagoon, green, light blue, lilac, lime, orange, pink, white, yellow



Variants

With side elastics (all colours) or with adjustable ties (light blue only)



Packaging

Pack of 50 masks with 'open/close' system
Box of 20 packs



CE Ref

Monoart® Face Mask Protection 3

Safe

Class I medical devices, designed to prevent airborne transmission of infectious agents, compliant with EN 14683 and ASTM F2100 standards.

Comfortable

Inner layer of hydrophilic Spunbond non-woven fabric with soft rounded elastics or adjustable ties and built-in nose bridge, combining comfort and moisture absorption in a Class I medical device

Practical

The "open/close" system ensures easy, hygienic access.

Face masks

Protection

Floral

Floral disposable medical face masks, fluid-resistant (21.3 kPa), protect against aerosols and splashes while preventing airborne microorganism transmission. High breathability (ΔP 39.2 Pa/cm²) exceeds EN 14683 standards. Features a 99.7% efficient antibacterial filter, soft rounded elastics, and a double-track nose bridge for comfort and stability.

Four-layer construction:

- Outer: double layer of coloured, water-repellent floral Spunbond fabric
- Middle: Meltblown non-woven filter
- Inner: hydrophilic Spunbond fabric for sweat and moisture absorption

Protective

Double outer layer of water-repellent Spunbond fabric for aerosol and splash protection, with an inner layer featuring particle and bacterial filters.

Fashion

Available in the original floral design and five vibrant colours, delivering a distinctive and easily recognizable appearance.

Comfortable

Inner layer of hydrophilic Spunbond non-woven fabric with soft rounded elastics and a double-track extended nose bridge, ensuring stable fit and enhanced comfort.



Classification

Medical Device Class I
EN 14683: type IIR
ASTM F2100: level 3



Colours

Floral burgundy, floral light blue, Floral lilac, floral lime, floral pink



Packaging

Pack of 50 masks with 'open/close' system
Box of 20 packs



CE Ref

Monoart® Face Mask
Protection 3 Floral

Face masks

Protection 4

Protection 4 disposable medical face masks, fluid-resistant (21.3 kPa), designed for aerosol and splash protection. Provide effective airborne microorganism filtration with high breathability (ΔP 44.3 Pa/cm²), exceeding EN 14683 standards. Equipped with a 99.7% efficient particulate and antibacterial filter, soft flat welded elastics, extended nose bridge, and Perfect Fit system for improved stability and comfort.

Four-layer structure:

- Outer: double layer of coloured water-repellent Spunbond fabric for enhanced protection
- Middle: Meltblown non-woven filter
- Inner: hydrophilic Spunbond fabric for moisture and sweat absorption



Classification

Medical Device Class I
EN 14683: type IIR
ASTM F2100: level 3



Colours

Black, green, light blue, lilac, orange



Packaging

Pack of 50 masks with 'open/close' system
Box of 20 packs



CE Ref

Monoart® Face Mask
Protection 4

Face masks FFP2 NR Protection

Monoart® FFP2 NR Protection face masks combine high protection, comfort, and style, available in six colours. Made in Italy with attention to detail, including colour-matched elastics and embossed branding, they offer technical excellence. These Category 3 PPE masks exceed EN 149:2001+A1:2009 standards with over 99% filtration and superior breathability, providing reliable respiratory protection against dust, non-volatile liquid particles, mists, and fumes.



Classification

Personal Protective Equipment Category 3, compliant with EU Regulation 2016/425



Colours

Black, blue, lilac, lime, orange, white



Packaging

Open/close pack with 10 individually wrapped face masks
Box of 20 packs



CE Ref

MONOART® FFP2 NR PROTECTION

Protective

Four-layer structure with Spunbond outer and inner layers and dual Meltblown filters, ensuring over 99% filtration efficiency.

Coloured

Available in six bright colours with tone-on-tone elastics for maximum protection without sacrificing style.

Comfortable

Soft, flat elastic bands welded to the outer part of the mask for maximum comfort and adherence to the user's face, all-Italian product quality.

Disposable latex gloves

Monoart® disposable latex gloves provide safe, hygienic hand protection with high comfort and durability. Featuring a polymer-coated inner layer and beaded cuff for easy donning, they ensure excellent tactile sensitivity and a secure grip in all conditions thanks to a micro-textured surface. Available in four colours for both performance and style.



Classification

Personal Protective Equipment Category 3, compliant with EU Regulation 2016/425



Colours

Blue, lilac, lime, pink



Size

XS, S, M, L



Packaging

Pack of 100 pieces
Box with 10 packs



CE Ref

MONOART® Latex Gloves

Safe

Provide a high-performance protective barrier, free from SVHC substances, endocrine disruptors, and dust particles commonly associated with latex-related respiratory reactions.

Certified

Manufactured under strict quality control and certified according to EU Regulation 2016/425, complying with EN ISO 374-1, EN ISO 374-4, and EN ISO 374-5 standards.

Resilient

Resistant to water, acids, and alkalis, providing effective protection against microorganisms and most caustic agents and detergents.

Disposable nitrile gloves

Monoart® disposable nitrile gloves are designed for latex-allergic users, providing a reliable protective barrier against contamination and infection. They deliver high mechanical strength and chemical resistance, suitable for procedures involving puncture risk, cuts, or aggressive substances. The micro-roughened surface ensures a secure grip, while maintaining excellent tactile sensitivity for precise handling.

High quality

Latex-free formulation, suitable for individuals with latex protein allergies, made from carefully selected materials to ensure reliable performance and quality.

Protective

Highly elastic and ergonomically designed, offering performance close to latex gloves. Suitable for use in situations with risk of puncture, chemical exposure, or viral contamination.

Resilient

Nitrile gloves are Category III PPE, offering enhanced mechanical strength and chemical resistance compared to standard latex and vinyl gloves.



Classification

Personal Protective Equipment Category 3, compliant with EU Regulation 2016/425



Colour

Black



Size

XS, S, M, L



Packaging

Pack of 100 pieces
Box with 10 packs



CE Ref

MONOART® Nitrile Gloves

Bandana

Cotton bandanas designed for comfort and easy wear, featuring adjustable ties for all sizes. Made from high-quality cotton with durable colour fastness, they are autoclavable at 121°C for reuse.

Product mix:

- Monoart® bandana set: 10 pieces in different colours (black, blue, light blue, light green, lilac, lime, orange, pink, white, yellow).



Classification

Personal Protective Equipment Category 1, compliant with EU Regulation 2016/425



Colours

Black, blue, green, light blue, light green, lilac, lime, orange, pink, white, yellow



Size

One size



Packaging

Individual pack
Box of 10 single-colour bags



CE Ref

Monoart® BANDANA

Comfortable

Made from soft cotton, designed for optimal fit and comfort even during extended use.

Autoclavable

Autoclavable at 121°C, made from high-quality cotton for reliable thermal resistance.

Universal

Fitted with adjustable ties for universal sizing, ensuring a secure and comfortable fit that conforms to the head shape.

Monoart® Infection Control Kit



Comprehensive infection-control kit designed to provide a reliable protective barrier for operators.

Contents include:

- 2 long gowns with zip closure
- 50 disposable polyethylene aprons
- 2 autoclavable cotton bandanas
- 50 Protection 3 face masks
- 2 pairs of Ultra Light protective glasses
- 2 pairs of cut-resistant gloves

Complete

The only all-in-one kit providing complete protection for operators throughout every stage of the sterilization process.

Certified

Kit of certified Medical Devices and PPE, crafted from high-quality materials and designed for professional use.

Reliable

Designed to reduce cross-contamination, this kit uses only certified, safe materials, backed by Monoart® quality.



Classification

Gown and gloves:
Personal Protective Equipment Category 3

Face Masks:
Medical Devices Class II
EN 14683: type II
ASTM F2100: level 1

Bandanas:
Personal Protective Equipment Category 1

Glasses:
Personal Protective Equipment Category 2



Packaging

1 box containing 1 kit



Gown

Full-length gown with zip closure, mandarin collar, and elastic cuffs for secure fit.



Apron

Full-length gown with zip closure, mandarin collar, and elastic cuffs for secure fit.



Cut proof gloves

Latex-free nitrile gloves, contoured and non-slip, offer comfort, sensitivity, and chemical resistance. Autoclavable, sizes S and M.



Classification

Personal Protective Equipment Category 3



Colour

White



Packaging

Carton of 25 pieces



Colour

White



Size

80 x 110 cm (l x h)



Packaging

Pack of 50 pieces
Box of 9 packs



Classification

Personal Protective Equipment Category 3



Colour

Violet



Packaging

Individual pack
Box of 9 packs

Protective glasses

Euronda Monoart® protective eyewear and visors are designed for operators and patients, providing reliable safety, ergonomic comfort, and modern design for clinical use.



Classification

Personal Protective Equipment Category 2, compliant with EU Regulation 2016/425



Variants

Available in different models and colours, with transparent or orange lens



Packaging

Individual pack
Box of 10 pieces (except Contemporary sold individually)



Contemporary



Cube
Cube Orange



Stretch
Stretch Orange



Evolution
Evolution Orange



FitUp



Light
Light Orange



Zero



Ultra Light



BabyOrange



Operator visor



Small Visor

Dental chair cover kit

Disposable dental chair cover set made from soft non-woven fabric (28 g/m²) for patient protection and hygiene. Includes headrest, backrest, and seat covers with elastic edges and ties for universal compatibility with all dental chairs.



Universal

Includes headrest cover, backrest cover and seat cover. Adaptable to all dental chair models thanks to the elastic edges and ties.

Comfortable

Made of soft, hypo-allergenic and breathable polypropylene non-woven fabric.

Hygienic

Protects the dental unit on a daily basis for patient safety.



Colours

Blue, green, light blue, lilac, lime, orange, white, yellow



Packaging

Pack of 1 kit consisting of 3 pieces
Box of 25 single-colour packs

Headrest cover

Headrest cover made of cellulose on the outside and polyethylene on the inside. Complete with practical dispenser with 'open/close' system.



Resilient

Maximum resistance and impermeability to protect the most sensitive area of the dental unit.

Universal

Its size and pouch design make it usable on various dental chair models for simple and immediate application.

Functional

The practical dispenser protects the product while ensuring hygiene and patient safety.



Colours

Blue, green, white



Size

28 x 30 cm



Packaging

Pack of 250 pieces
Box with 8 packs

Tray Paper

Steam- and gas-permeable crepe paper, free of chlorine, dioxin and ash. Developed by Euronda research to support the operation and sterilization process, its paper does not degrade or crumble on contact with steam, making it the ideal liner for the tray.

Product assortment available

- Mix Monoart® Tray Paper: contains 6 packs of Tray Paper (light blue, lilac, lime, orange, pink, yellow).



Colours

Green, light blue, lilac, lime, orange, pink, white, yellow



Size

28 x 18 cm



Packaging

Pack of 250 pieces
Box with 10 packs

Practical

Tear-resistant, highly absorbent crepe paper sheets. The ideal liner for the bottom of the tray.

Resilient

Steam- and gas-permeable and resistant to moisture, they do not fray or degrade during sterilization, avoiding residues that could clog the autoclave drain.

Safe

Free of chlorine, dioxin and ash, they minimize water condensation.

Decontamination mat

Monoart® decontamination mat in low-density polythene with 30 numbered and branded adhesive sheets, to prevent the spread of bacteria in the surgery. Thanks to the special formula of the bactericidal gelatine present on the upper surface, it removes and deactivates the bacteria present on the soles of shoes. Product made in Italy.



Classification

Medical Device Class I



Colour

Blue



Size

45 x 115 cm



Packaging

Pack of 4 pieces, each one with 30 removable adhesive sheets



CE Ref

Euronda MONOART - Tappeto Decontaminante Multistrato, Euronda MONOART - Multilayer Decontaminating Mat

Practical

The numbered sheets optimize the handling of the decontamination mat by assisting the operator in changing the different layers.

Hygienic

An additional layer of protection that increases the level of safety in operating theatres and surgeries, supporting the creation and maintenance of the sterile field.

Functional

Smooth surface without imperfections, for maximum effectiveness in removing dirt. Fixing band at the bottom for perfect grip to the floor.

Barrier film

Disposable film for the protection of areas and surfaces, especially those difficult to reach, clean or disinfect. Reduces the risk of cross-contamination caused by aerosols and sprays. Thanks to its lightweight adhesive, it adheres easily to all surfaces without leaving residue, while the non-adhesive border on both sides makes placement and removal quick and easy.



Colour

Transparent blue



Size

Pack 18 x 12 x 16.5 cm (l x d x h)
Sheet 10 x 15 cm



Packaging

Single pack containing one roll of 1200 tear-off sheets with dispenser

Safe

An all-purpose infection-control barrier designed to isolate all instruments and surfaces quickly and easily.

Protective

Quick to apply and remove, it leaves no residue and is suitable for use on all surfaces and objects.

Practical

Pack contains a roll with 1200 tear-off sheets and is designed to ensure the highest level of hygiene and ease of use.

Turbine cover

Disposable sleeve in transparent polyethylene film, which helps prevent the risk of cross-contamination. Product made in Italy.

Available in two different sizes:

- 4.5 cm diameter, ideal for turbine cords, suction tubes and handpieces;
- 5.5 cm diameter for micromotor cords, air/water syringes and connecting tubes.



Colour

Transparent



Size

Pack 11 x 6 x 11 cm (l x d x h)
Sheets:
4.5 x 45 cm (Ø x l)
5.5 x 45 cm (Ø x l)



Packaging

Single pack containing one roll of 250 tear-off sheets

Functional

Guarantees and simplifies the protection of devices and components used in close contact with the patient.

Practical

Pack containing one roll of 250 tear-off sheets and a dispenser, designed to ensure maximum hygiene and ease of use.

Tailor-made

Available in two sizes with a diameter of 4.5 and 5.5 cm, designed to optimally protect various suction systems and micromotors.

Dappen dishes

Practical paste holders for cleaning and polishing teeth as prophylaxis, they can also be used for holding composite and etching fluid and can be used to mix dental cements. Available in assorted shades so you can personalize your dental practice with a touch of colour and originality.

Product mix available

- Mix Monoart® Dappen dishes: contains 4 packs of dappen dishes (50 pieces) per colour (blue, fuchsia, lilac, lime, orange).



Colours

Blue, fuchsia, lilac, lime, orange



Size

3.3 x 2.4 cm (Ø x h)



Packaging

Pack of 50 pieces
Box with 20 packs

Functional

Paste containers made of polypropylene, ideal as composite holders, for etching fluids and for mixing cements.

Ergonomics


Equipped with slots, enabling you to attach them to your tray, and the grooved surface ensures a secure grip and makes them easy to handle.

Coloured

Available in five assorted shades so you can personalize your dental practice with a touch of colour and originality. Try out the five colours with the Dappen dishes Mix for a different shade each day!

Bur blocks

Bur blocks for dentists with perforated base and transparent lid, suitable for storing diamond burs, contra-angle and handpiece. The base is autoclavable. Available in different formats.



Colour

Light grey



Size

Bur holders blocks with autoclavable base for:
18FG, 18HS or mixed 18 FG RA

Bur holders blocks with autoclavable base for:
36FG, 36HS or mixed 36 FG HS

Bur holders blocks with autoclavable base for:
72FG, 72HS or mixed 72 FG HS

Simple

Essential design, compact dimensions and polypropylene base autoclavable to a maximum of 130°C.

Versatile

Available in different formats, 18, 36, 72 bur blocks with FG, HS or mixed holes, to meet all different requirements.

Resilient

Made of high-quality materials, they are durable and long-lasting to support the daily work of the dental practice for a long time.

Alginate mixing bowls

Bowls for mixing dental material before using on the patient. Available in different formats. Product made in Italy.



Functional

Non-toxic and latex-free PVC cups, ideal for mixing plaster and alginates.

Versatile

Available in four different formats, they adapt to the needs of every dental practice.

Practical

The smooth surface makes cleaning the cups extremely easy and straightforward. The ergonomic shape makes it easy to use.



Colour

Green



Size

Bowl No. 1 (small)
11 x 7.5 cm (Øxh)

Bowl No. 2 (medium)
12.5 x 9 cm (Øxh)

Bowl No. 3 (large)
14 x 11.2 cm (Øxh)

Bowl No. 4 (maxi)
16.5 x 15.4 cm (Øxh)



Packaging

Single pack

Spray ice

Coolant spray ice in 200 cc can. Ideal for cooling impression material such as hydrocolloids and glass plates.



Classification

Medical Device Class I



Shelf life

3 years



Packaging

Box of 12 cans



CE Ref

Cold Spray Classic

Functional

Hypothermizing spray consisting of a mixture of hydrocarbons. Reaches a maximum temperature of -45°C.

Practical

Easy to use thanks to the precision nozzle.

Safe

Non-toxic product with a delicate mint fragrance.

CDS301 stool

Chromed steel frame with 5-spoke base, equipped with five self-propelled castors. The convenient mechanism is integrated into the underside of the seat and can be operated 360° using both hands, so that searching for levers and buttons during procedures is no longer necessary. The wraparound, anatomical backrest is adjustable in height, depth and position. It allows you to maintain a relaxed working position, reducing the load on your back, pelvis and legs. Smooth surfaces for maximum hygiene and cleanliness. Product made in Italy.



Colours

Numerous colours available



Models

Thermoformed or sewn



Size

Round seat
40 x 10 cm (Øxh)
Adjustable height
51 – 69 cm



Weight

15 kg



Coral stool

Operator stool with recyclable, fire-retardant, scratch- and impact-resistant extruded polyamide frame, equipped with five self-locking castors. The operating mechanism is covered by a special cap. It is integrated into the underside of the seat and can be operated 360° using both hands, so that searching for levers and buttons during procedures is no longer necessary. The wraparound, anatomical backrest is adjustable to three different positions. The seat is made of non-deformable polyurethane and covered with sky. Product made in Italy.



Colours

Numerous colours available



Models

Thermoformed or sewn



Size

Round seat
39 x 9 cm (Øxh)
Adjustable height
52 – 70 cm



Weight

13 kg



Onyx stool

Aluminium frame with five spokes. Adjustable height thanks to the lever located under the seat, which can be operated 360°. Round seat and wraparound, anatomical backrest, adjustable to three positions (right, centre and left) and by depth. Product made in Italy.



Colours

Numerous colours available



Models

Thermoformed or sewn



Size

Round seat
39 x 9 cm (Øxh)
Base diameter 58 cm
Adjustable height
50 - 65 cm



Weight

11 kg







| Parque Empresarial do Padrão, Lote11, 5460-343 Boticas, Portugal |
| +351 276 418 287 | info@sterifast.com | www.sterifast.com |